



Track & Trace, Authentication System





- FO-Security was founded in 1930
- Guilloche machine "Fotorotar" sold worldwide
- car registration cards, ID cards, vouchers, public transport tickets, stock and bond certificates, security labels, serialized packaging
- EN-ISO 9001, EN-ISO 14000, EN-ISO 27001, OHSAS 18001, FSC IMO-COC-028635
- member of Intergraf
- FO Group: 25 mio CHF annual sales, 120 employees



FO-Security, Ein Geschäftsbereich der FO Print & Media AG | Gewerbestrasse 18 | 8132 Egg/ZH | Tel +41 44 986 35 40 | Fax +41 44 986 35 36 | E-Mail info@fo-security.ch | www.fo-security.ch



FO-Security in partnership with JURA

- JURA is a leading software developer for security pre-press
- CWA 15374:2005 certificate, based on ISO 9001:2000, by Intergraf (International Confederation for Printing and Allied Industries) and VPGI
- graphic software for security document design
- Direct Laser Engraving (DLE®) for intaglio plates
- CTP output devices developed in co-operation with Kodak
- personalization solutions for travel and personal documents
- covert security features for banknotes, security documents

FO-Security, Ein Geschäftsbereich der FO Print & Media AG | Gewerbestrasse 18 | 8132 Egg/ZH | Tel +41 44 986 35 40 | Fax +41 44 986 35 36 | E-Mail info@fo-security.ch | www.fo-security.ch









- Uses a 2D code, works with multiple types
- Does not require special equipment, can be printed with regular inkjet printers
- Tracking and authentication happen as one step, at the same time
- Uses a smartphone with a reader application
- Automatically logs all events and results, which are available live for customers



Step 1: A serialised 2D code is generated , including product code, serial no, expiry date and batch number







Step 2: The 2D code is input into IQ-R encoding system.



- Encodes a second layer of information.
- The secondary information is hidden from public reader applications.
- The specific hidden information is entered into a database.



Step 3: The new IQ-R code contains both the public primary and the hidden secondary information and is ready for printing



The restulting IQ-R is indistinguishable from the original by the bare eye or any reading application.



Step 4: Authentication requires nothing more than a smartphone application


- The application can identify on spot if the product (2D code) is genuine or copy
- The application can contact a server to download the code-specific hidden information and to compare it to the readout





Step 5: Based on the analysis, the application provides feedback whether the IQ-R code is authentic









Genuine Code

Improper read out conditions, retry required Forged code




F SECURITY

















The process of tracking, tracing and authentication




- Uploads data when authenticating the product
- Downloads data for authenticity check

Benefits:

- Can discover counterfeited medicine before buying them
- Knows that forging happened before reaching this step
- Provides feedback to the pharmacy
- Provides feedback to the manufacturer via the live information in the database
- Patient can be targeted with extra services via the application







- Contains product code, serial no, expiry date, batch no
- Easy to implement; no specialised tools for manufacturing
- Does not require specialised devices for tracking and authenticating
- Provides protection against forgery
- Provides protection against photocopying
- Provides live tracking and authentication information
- Can be used as a channel to reach the patients with other types of information, such as:
 - R&D Doctor Patient relationship management
 - Package insert





